

original claim  
staff. a c  
is a small  
May 6

[illegible] $\beta$ 

/

The tests leading to the establishment of the primer sequences for the BRCA1 and hMLH1 of the present invention were conducted with the TDGS design prepared with the computer programming and equipment described in PCT/IB97/00976, published on or about February 14, 1998.

### Objects of Invention and Summary

The objects of the invention are to provide novel theoretically and empirically (experimentally) derived TDGS patterns for hMLH1 and BRCA1 genes which may be used by testers to test for gene sequence variation and/or mutations.

### Drawings

Figs. 1A and 1B show the computer-aided design TDGS patterns obtained for the hMLH1 and BRCA1 (theoretical-left hand side; empirical or experimental--right hand side).

In the theoretical vs. empirical patterns of the MLH1 and BRCA1 genes, for all four genes, one or more exons were designed in overlapping fragments, in which case the fragment name is exon.1, exon..2, etc. Exons 8 and 15 of hMLH1 contain polymorphisms, which can be distinguished from disease-causing heterozygous mutations on the basis of a unique four-spot pattern (18).

### Description Of The Invention In Preferred Forms

*The MLH1 DNA mismatch repair gene.* The design for *MLH1* took 30 minutes (excluding exon indication). Fig. 1A shows the theoretical and the empirical TDGS pattern for the *MLH1* gene. Because exons 11 and 12 had to be subdivided into overlapping fragments, two multiplex groups are currently being used, with the long PCR carried out as a four-plex PCR. Like many other genes, exon 1 of *MLH1* is GC-rich and, hence, was found to melt at a much higher % UF compared to most of the other fragments. Thus far, a total of 41 coded samples with previously identified mutations have been analyzed in a blinded fashion with 100% concordance (30).

*The breast and ovarian cancer susceptibility gene BRCA1.* The tumor suppressor gene *BRCA1* contains 24 exons, of which exon 11 contains approximately 60 % of the coding region. Fig. 1B shows the theoretical and empirical 2-D pattern for *BRCA1*. Of all 2-D designs discussed, this was the most difficult (total design time was 2 h), the main reason being the need to make overlapping fragments for the 3.4 kb exon 11. Pre-amplification was accomplished by one 7-plex long PCR. Using the long PCR amplicons as template, all 24 exons were amplified in a total of 37 fragments distributed over 5 multiplex groups. The overlap and sometimes short distances from fragment to fragment necessitated the use of so many multiplex short PCR groups. The non-coding exons 1a, 1b and the non-coding part of exon 24 were excluded. Evaluation of this test design using a panel of coded samples with previously identified mutations is currently ongoing. Thus far, mutations and polymorphisms have been detected in exons 2, 8, 11, 16, 20 and 23.

### PCR Amplification

Primers were obtained from Genosys Biotechnologies, Inc. (The Woodlands, TX). For complete lists of all sequences, except *BRCA1*, see references 18, 29 and 30. Primer sequences for *BRCA1* will be published elsewhere but will be made available upon request. PCR amplification of gene sequences was carried out using the two-step protocol first described by Li and Vijg (22). Primers for long-distance PCR were designed based on published sequences (24-27) using Primer Designer 3, to amplify the entire gene-coding region for each of the 4 genes as a 1-plex (*TP53*), a 6-plex PCR (*RBI*), a 4-plex PCR (*MLH1*) or a 7-plex PCR (*BRCA1*). The LA PCR kit (Takara) was used for long PCR in a PTC-100 thermocycler (MJ Research). Multiplex short PCR was carried out using the long PCR products as template. Between 0.1 and 1.125  $\mu\text{M}$  of each primer was used in a 50  $\mu\text{l}$  reaction with 1  $\mu\text{l}$  of long PCR product in 20 mM Tris-HCl (pH 8.4), 50 mM KCl, 250  $\mu\text{M}$  of each dNTP and 5 % formamide. Two and a half units of Taq DNA polymerase (Life Technologies) were added after an initial denaturation at 94 °C for 60 s. Cycling conditions for multiplex short PCR and concentrations of  $\text{MgCl}_2$  varied among different genes and amplifications were carried out in a PTC-100 thermocycler (MJ Research).

### Two-dimensional DNA electrophoresis

For *RBI*, 5  $\mu\text{l}$  of multiplex short PCR was used per electrophoresis run. For *TP53*, *MLH1* and *BRCA1*, 5  $\mu\text{l}$  of each of the different multiplex groups were combined. One tenth of a

5

-5  
volume of loading buffer (0.25 % xylene cyanol, 0.25 % bromophenol blue, 15 % ficoll and 100 mM Na<sub>2</sub>EDTA) was added and the mixtures were loaded onto a 6.5 % (TP53) or 10 % (RBI, MLH1 and BRCA1) PAA non-denaturing size gel (acrylamide: bisacrylamide = 37.5:1) in 0.5 x TAE buffer. The samples were electrophoresed for 5.3 h at 150 V (RBI), 5 h at 120 V (TP53) or 7.5 h at 140 V (MLH1 and BRCA1) at 50 °C. After staining the gel with a mixture of equal amounts of SYBR-green I and II (Molecular Probes, Eugene, Oregon) for 20 min, the region containing all fragments of interest (usually between 100 and 600 bp) was cut out and loaded onto a denaturing gradient gel (DGGE). Gradients used were 0 to 50 % UF for RBI, 20 to 70 % UF for TP53, 25 to 70 % UF for MLH1 and 20 to 65% UF for BRCA1. The second orthogonal dimension was run for 12 h at 100V (RBI), 14 h at 120 V (TP53) or 16 h at 100 V (MLH1 and BRCA1). Spot patterns were visualized by SYBR-green staining using a FluorImager (Molecular Dynamics, Sunnyvale, California).

069050-00000000

6 The primer sequences for long and short PCR for the BRCA1 are as

follows:

A. Primer Pairs for Long Distance PCR

Exons 1-4

MLH1-4F GCG.GCT.AAG.CTA.CAG.CTG.AAG.GAA.GAA.CGT.GA

MLH1-4R GGC.GAG.ACA.GGA.TTA.CTC.TGA.GAC.CTA.GGC.CC

product size= 10.8kb

Exons 5-10

MLH5-10F

GCG.CCC.CTT.GGG.ATT.AGT.ATC.TAT.CTC.TCT.ACT.GG

MLH5-10R GCG.CTC.ATC.TCT.TTC.AAA.GAG.GAG.AGC.CTG

product size=10.5kb

Exons 11-13

MLH11-13F GCG.CTT.TTT.CTC.CCC.CTC.CCA.CTA.TCT.AAG.G

MLH11-13R GGG.TTA.GTA.AAG.GAA.GAG.GAG.CTT.GCC.C

product size=8.7kb

Exons 14-19

MLH14-19F GGT.GCT.TTG.GTC.AAT.GAA.GTG.GGG.TTG.GTA.G

MLH14-19R

GCG.CGC.GTA.TGT.TGG.TAC.ACT.TTG.TAT.ATC.ACA.C

product size=10.5kb

Underlined nucleotides represent nucleotides added to modify melting temperatures of the primers

B. Primer Pairs for Short PCR

Exon Clamp<sup>1</sup> Product Size Tm<sup>2</sup> Primer Sequence

12.1	40	184	44.53	TTT.TTT.TTT.TTT.TAA.TAC.A AAT.CTG.TAC.GAA.CCA.TCT
12.2	8	366	53.23	TGG.AAG.TAG.TGA.TAA.GGT TGT.ACT.TTT.CCC.AAA.AGG
	40			
13	40	272	49.06	ATC.TGC.ACT.TCC.TTT.TCT AAA.ACC.TTG.GCA.GTT.GAG
14	45	235	48.94	TAC.TTA.CCT.GTT.TTT.TGG GTA.GTA.GCT.CTG.CTT.GTT
	5			
15	40	179	29.97	CAG.CTT.TTC.CTT.AAA.GTC CAG.TTG.AAA.TGT.CAG.AAG
16		261	47.56	CTT.GCT.CCT.TCA.TGT.TCT.TG AGA.AGT.ATA.AGA.ATG.GCT.GTC
	40			
17	40	199	47.01	A'TT.ATT.TCT.TGT.TCC.CTT AAT.GCT.TAG.TAT.CTG.CCT
18	45	215	46.67	CCT.ATT.TTG.AGG.TAT.TGA.AT GCC.AGT.GTG.CAT.CAC.CA
19.1		282	43.43	TGT.TGG.GAT.GCA.AAC.AGG ATC.CCA.CAG.TGC.ATA.AAT
	40			

1 GC clamps:

50 clamp:

CGC.CCG.CCG.CCG.CCC.GCC.GCG.CCC.CGC.GCC.CGT.CCC.GCC.GC  
C.CCC.GCC.CG

45 GC clamp:

*clamps 27-32*

669050-EE90E60

CGC.CCG.CCG.CGC.CCC.GCG.CCC.GTC.CCG.CCG.CCC.CCG.CCC.GG  
C.CCG

40 clamp:

CGC.CCG.CCG.CGC.CCC.GCG.CCC.GGC.CCG.CCG.CCC.CCG.CCC.G

8 clamp:

CGT.CCC.GC

5 clamp:

GCG.CG

2 clamp:

CG

<sup>2</sup>T<sub>m</sub> is given in %UF

669090 44930260



### Primers for long-PCR BRCA1 (7-PLEX)

**Size:** 9.9 kb

**Size:** 9.7 kb

**Size:** 4.8 kb

**Size:** 9.0 kb

**Size:** 10.7 kb

**Size:** 7.2 kb

**Size:** 11.4 kb

[illegible]

## BC1EX11

Exon Frag Primers 5' -&gt; 3'

		size	Tm(%UF)
11.1	[GC 3]ACCTTGTTATTTTGTATATTT 22 [GC 13]TTGCTAAGCCAGGCTGTT 18	347	40.99
11.2	[GC 3]ATACTCATGCCAGCTCATT 20 [GC 12]AACGTCCAATACATCAGCTA 20	461	40.74
11.3	CATGCTCAGAGAATCCTAGA 20 [GC 3]CTGTGGCTCAGTAACAAATG 20	438	35.04
11.4	[GC 12]TCACTCCAAATCAGTAGAGA 20 [GC 3]TACTGCTGCTTATAGGTTCA 20	476	34.85
11.5	[GC 3]GAAAGCAGATTGCGCAGTTC 20 [GC 11]CTGACTGGCATTGCGTTGTA 20	468	33.66
11.6	[GC 3]GAATAGGCTGAGGAGGAAGT 20 [GC 13]CTCTTGGAAGGCTAGGATTG 20	410	40.51
11.7	[GC 3]ACAGCGATACTTTCCCAGAG 20 TGCCTTCCCTAGAGTGCTAA 20	345	36.45
11.8	TTGCCAAACTGAAAGATCTGT 20 [GC 3]GCTTTGAAACCTTGAATGTA 20	365	38.37
11.9	[GC 13]GTCGGGAAACAAGCATAGAA 20 [GC 4]TTGCCTCTGAACTGAGATGA 20	422	40.40
11.10	[GC 12]TAATATCACTGCAGGCTTTC 20 [GC 1]TTCCTCAAAGTTTTCCTCTA 20	292	35.93
11.11	[GC 1]TCCCATCAAGTCATTGTTA 20 TTCCAGGAAGACTTTGTTTA 20	390	33.06
11.12	[GC 12]TAATGAAGTGGGCTCCAGTA 20 [GC 1]CTTCCCATAGGCTGTTCTAA 20	309	33.22
11.13	[GC 1]GCAAGAATATGAAGAAGTAG 20 CAAATGTGTATGGGTGAAAG 20	305	37.43
11.14	[GC 1]AGACACCTGATGACCTGTTA 20 [GC 12]TCTCCTCTGTGTTCTTAGAC 20	378	43.03
11.15	CCTTTCACCCATACACATTT 20 [GC 8]GACTGATGCCTCATTTGTTT 20	460	39.33
11.16	[GC 3]CTCAGGAACATCACCTTAGT 20 [GC 16]ATAAATAGACTGGGCCACAC 20	356	44.00

669690 E E 5050

## All exons excluding exon 11

BRCAONE

Exon Frag Primers 5' -&gt; 3'

size Tm(%UF)

2	1	[GC 1]TATATATGTTTTCTAATGTGT 22 [GC 12]TCCCAAATTAATACACTCTT 20	203 34.64
3	1	[GC 12]GAGCCTCATTATTTTCT 18 [GC 4]ATTTTTCGTTCTCACTTA 18	269 37.22
5	1	[GC 4]TATTGCGCTTTTGAGTAT 18 [GC 12]TCTGATGAATGGTTTTAT 18	305 26.69
6	1	[GC 8]ACTTGCTGAGTGTGTTTC 18 GCACTTGAGTTGCATTCT 18	213 35.52
7	1	[GC 3]TACATTTTTCTCTAACTGC 19 GAAGAAAACAAATGGTTTT 19	250 32.67
8	1	GGAGGAAAAGCACAGAAC 18 [GC 3]CCAGCAATTATTATTAAATACTT 23	248 40.51
9	1	[GC 3]CAGTAGATGCTCAGTAA 18 AATACCAGCTTCATAGAC 18	242 24.26
10	1	[GC 4]CTGCATACATGTAAC TAG 18 CTACCCACTCTCTTTTCA 18	229 38.30
12	1	[GC 4]AGTTGCAGCGTTTATAGT 18 [GC 13]CAGCAAACCTAAGAATGT 18	289 48.54
13	1	[GC 4]GCTTCTCAAAGTATTTC A 18 AGTGTTTGGCCAACAATA 18	293 45.18
14	1	[GC 4]CCAATTTGTGTATCATAG 18 [GC 13]AGTGTATAAATGCCTGTA 18	417 30.78
15	1	[GC 1]TGGTTTTCTCCTTCCATTTA 20 [GC 16]TGTTCCAATACAGCAGATGA 20	303 46.07
16	1	[GC 13]CGTTGTGTAAATTAACTTC 20 [GC 1]AGTCATTAGGGAGATACATA 20	427 47.49
17	1	[GC 4]TGTGCTAGAGGTA ACTCA 18 [GC 11]CTCATGTGGTTTTATGCA 18	242 32.51
18	1	[GC 12]TTTCAACTTCTAATCCTTT 19 [GC 4]GGAGAAATAGTATTATACT 19	194 36.32
19	1	GTCTTCTGCTGTATGTA 18 [GC 4]CTGAATGAATATCTCTGG 18	178 32.32
20	1	[GC 4]CTCTTTCTCTTATCCTGAT 19 TGGTGGGGTGAGATTTTT 18	219 46.40
21	1	[GC 8]ATTCCCTGTCCCTCTCT 18	172 49.95

CTGGAACTCTGGGGTTCT 18

2 1 [GC 4]TGATTTTACATCTAAATGTC 20  
[GC 13]AGGAGAGAATATTGTGTC 18

209 47.71

3 1 [GC 12]TAGTCCTACTTTGACACT 18  
[GC 4]AAATATTTAAATGTGCCAA 20

275 49.47

4 1 [GC 13]AATCTCTGCTTGTGTTCTCT 20  
[GC 18]ATTTAGTAGCCAGGACAGTA 20

325 59.79

005050-EE9DE60